

Cancer Antigen 27.29

Effective Date: May 12, 2017

All Cancer 27.29 antigen testing for Beaumont Laboratory are performed at Warde Medical Laboratory. Effective May 12, 2017 Warde Medical Laboratory will be unable to perform dilutions on samples with CA27.29 results greater than 450 U/mL as this dilution protocol has been discontinued by the manufacturer of the assay (Siemens). Since all major referral labs use the Siemens platform for testing, no alternate method is available at this time. Warde Medical Laboratory will resume dilutions for the assay as soon as the manufacturer validates new diluent. In the meantime, all results greater than 450 U/mL will be resultated as such, without any further quantitation.

If you have questions, please contact Client Services (1-800-551-0488, option 5).

Laboratory Test Directory: <http://beaumontlaboratory.com/test-lab-directory>.

Date submitted: 4/25/2017

Submitted by: Yvonne Posey, M.D., Associate Medical Director, Automated Chemistry & Special Testing