

VITAMIN E

Effective Date: September 5, 2018

Effective September 5, 2018, specimens submitted for Vitamin E ***must be light protected.***

In addition, ***room temperature stability is no longer acceptable.***

Temperature	Old Stability	New Stability
Room Temperature	14 days	Unacceptable
Refrigerated	14 days	14 days
Frozen	14 days	14 days

If you have questions, please contact Client Services (1-800-551-0488, option 5).

For additional information, access the Laboratory Test Directory at:

<http://beaumontlaboratory.com/test-lab-directory>

Date submitted: September 14, 2018

Submitted By: Yvonne Posey, MD, Medical Director, Send Outs, Royal Oak

Beaumont Laboratory
Customer Service
1-800-551-0488