

CHANGE IN SPECIMEN REQUIREMENT FOR COMPREHENSIVE DRUG ANALYSIS

Effective Date: 12/04/2018

The laboratory will no longer require a blood specimen for the Comprehensive Drug Analysis with Confirmation testing, as all testing will be performed on urine specimens.

Acetaminophen and salicylate will be screened for on the urine specimen. If the initial screen is positive, the result will be accompanied by a comment that states:

- Acetaminophen detected in urine. If acute overdose is suspected, serum level is suggested.
- Salicylate detected in urine. If acute overdose is suspected, serum level is suggested.

Test Name	Comprehensive Drug Analysis with Confirmation
SOFT/EPIC Codes	CMPDA/LAB6854
Specimen Collection	Random urine
Specimen Volume	50 mL urine (30 mL minimum, 20 mL for pediatrics)
Specimen Transport	Refrigerated (2-8 °C)
Performed	Routine: Results available within 24 hours. STAT. Results available within 3 hours
Test Methodology	Immunoassay, Gas Chromatography/Mass Spectrometer (GC/MS), Liquid Chromatography w/ Tandem Mass Spectrometry (LC/MS/MS)
CPT Codes	80307x8, 80304, 80329, 80345, 80347, 80349, 80353, 80358, 80364, 80375, 83992 (alternate: G0479, G0480, G0482).
Test Name	Comprehensive Drug Analysis with Confirmation
SOFT/EPIC Codes	CMPDA/LAB6854

If you have questions, please contact Client Services (1-800-551-0488, option 5).

Laboratory Test Directory: <http://beaumontlaboratory.com/test-lab-directory>.

Date submitted: 11/20/2018
Submitted by: Michael Smith, PhD, Technical Director, Toxicology and Therapeutic Drug Monitoring Laboratory.
Timothy Kennedy, MD, Associate Medical Director, Automated Chemistry

Beaumont Laboratory
Customer Service
1-800-551-0488