

Beaumont Laboratory's tube chart shows the most commonly ordered tests with their associated specimen tube. Please confirm proper preparation and most updated tube type in our Laboratory Test Directory located at beaumontlaboratory.com or call our client service department at **800-551-0488**.

SERUM SEPARATOR TUBE (SST)	Acute Hep Pl.	BUN	FSH	HIV 1/2 Ab	LH	PTH ^{•◇}	TSH
	AST	Creatinine	HBcAB, Total	Homocysteine	Magnesium	PSA	Uric Acid
	ALT	CRP-HS	HBsAg	HSV I & 2, IgG	Phosphorus	Renal Func. Pl.	Vit B-12
	Amylase	Comp. Met. Pl.	HBsAb	Iron/TIBC	Potassium	Rheum. Factor	Vit D, 25 Hydrox
	Alk. Phos.	Electrolytes Pl.	HCV Ab	LD	Progesterone	RPR	
	Basic Met Pl.	Ferritin	hCG Beta Quant	Lipase	Protein Electro	T4, Free	
	Bilirubin	Folic Acid	Hepatic Func. Pl.	Lipid Pl.	Protein, Total	Testosterone	

LAVENDAR TUBE <i>K₂ EDTA</i>	BNP	Hemoglobin and	Path Consult
	CBC/Plt	Hematocrit	Retic Count
	CBC/Diff & Plt	Hemoglobinopathy	Sed Rate
	Factor V Leiden	Evaluation	Vit B-6 ^{•*}
	Hgb A1C	HIV Quant	

LT. BLUE <i>3.2% Sodium Citrate</i>	Anti-Thrombin III [•]	Factor V Activity [•]	Inhibitor Screen [•]	PTINR
	D-Dimer [•]	Hexagonal Phase	Protein C Activity [•]	PTT
	DRWT [•]	Phospholipid [•]	Protein S Activity [•]	

RED TUBE <i>Plain</i>	Carbamazepine	Theophylline
	Digoxin	Tobramycin
	Phenobarbital	Valproic Acid
	Phenytoin	Vancomycin
	Phenytoin, Free	

PINK TUBE <i>EDTA</i>	ABO & RH Type
	Ab Screen, RBC
	ABO
	RH Type

GREY TUBE <i>Potassium Oxalate/ Sodium Flouride</i>	Glucose Tolerance Test	Glucose Tolerance 2hr, Pregnancy
	Glucose Pregnancy Screen, 1hr	Glucose Tolerance 3hr, Pregnancy

GREEN TUBE <i>Sodium Herapin</i>	Vit B-1 (Thiamin) [*]
--	--------------------------------

LT. GREEN TUBE <i>Lithium Herapin Plasma Gel. Sep.</i>	Troponin I
--	------------

ROYAL BLUE TUBE <i>EDTA</i>	Lead Level, Whole Blood	Mercury
	Heavy Metals Screen	Zinc RBC
	Magnesium, RBC	

ROYAL BLUE TUBE <i>Plain</i>	Aluminum	Nickel
	Copper	Zinc Level
	Manganese	

• frozen pour over - special preparation * protect from light ◇ PTH requires 1 frozen pour over and 1 SST refrigerated